

Elabscience®

Annexin V-APC/7-AAD 细胞凋亡检测试剂盒

· 助力更好检测细胞凋亡 ·




APC与7-AAD之间荧光干扰小，数据分析复杂性大大降低



7-AAD同PI有着相似的荧光特性，但发射波谱比PI窄，对其他检测通道的干扰更小



浓缩型/即用型两种Binding Buffer可供选择，储存/实验操作更方便



Jurkat细胞用5 μ m喜树碱处理4 h后，用Annexin V-APC / 7-AAD 细胞凋亡检测试剂盒(E-CK-A218)检测结果

更多7-AAD凋亡检测试剂盒

货号	名称	规格
E-CK-A212	Annexin V-FITC/7-AAD Apoptosis Kit	20 / 50 / 100 / 200 Assays
E-CK-A216	Annexin V-PE/7-AAD Apoptosis Kit	20 / 50 / 100 / 200 Assays
E-CK-A218	Annexin V-APC/7-AAD Apoptosis Kit	20 / 50 / 100 / 200 Assays
E-CK-A218B	Annexin V-APC/7-AAD Apoptosis Kit	20 / 50 / 100 / 200 Assays

*注：E-CK-A218B产品内缓冲液为1 X Annexin V Binding Buffer

Elabscience Biotechnology Co.,Ltd.

Elabscience®提供全面准确的Annexin V细胞凋亡检测试剂盒



原料自研，品质溯源，价格更实惠



操作简单，15~20 min即可完成染色



性价比高，100 Assays可检测200次



15种荧光素搭配3种核染料，35种荧光搭配试剂盒，满足您的个性化需求



Elabscience®2022-2023年Annexin V高分文献展示（部分）

产品	文献信息
E-CK-A211 (Annexin V-FITC/PI)	MARCKSL1-2 reverses docetaxel-resistance of lung adenocarcinoma cells by recruiting SUZ12 to suppress HDAC1 and elevate miR-200b (<i>Molecular Cancer</i> , IF=41.44)
E-CK-A254 (Annexin V-Elab Fluor® 647/DAPI)	Cell-in-Cell-Mediated Entosis Reveals A Progressive Mechanism in Pancreatic Cancer (<i>Gastroenterology</i> , IF=29.4)
E-CK-A217 (Annexin V-APC/PI)	Cell membrane-camouflaged bufalin targets NOD2 and overcomes multidrug resistance in pancreatic cancer (<i>Drug Resistance Updates</i> , IF=24.3)
E-CK-A217 (Annexin V-APC/PI)	Pyroptosis Remodeling Tumor Microenvironment to Enhance Pancreatic Cancer Immunotherapy Driven by Membrane Anchoring Photosensitizer (<i>Advanced Science</i> , IF=17.521)
E-CK-A211 (Annexin V-FITC/PI)	Lysine methylation promotes NFAT5 activation and determines temozolomide efficacy in glioblastoma (<i>Nature Communications</i> , IF=16.6)



明星产品推荐

Annexin V-FITC/PI 细胞凋亡检测试剂盒

自2019年7月份第一篇文献引用，截止到2024年9月，该明星产品已发表文献**380+**篇。有一大批高质量论文发表在*Cells*, *Molecular Cancer*, *Journal*, *Nature Communications*, *Biomaterials*等国际知名期刊，最高影响因子**IF=41.44**。

货号	名称	规格
E-CK-A211	Annexin V-FITC/PI Apoptosis Kit	20 / 50 / 100 / 200 Assays
E-CK-A211B	Annexin V-FITC/PI Apoptosis Kit	20 / 50 / 100 / 200 Assays

*注：E-CK-A211B产品内缓冲液为1 X Annexin V Binding Buffer

☎ 免费技术热线：400-9627-365

武汉伊莱瑞特生物科技股份有限公司

Elabscience Biotechnology Co., Ltd.

☎ 电话：400-999-2100

🌐 官网：www.elabscience.cn

✉ 邮箱：marketing@elabscience.cn

📍 地址：武汉市高新大道858号光谷生物医药产业园二期B18栋



企业官网



微信公众号

EFL75-2410v2