


ELISA Kits

for Digestive System Diseases Research



//

慢性胃炎/肠炎

药物性肝损伤

肝硬化

急/慢性腹泻

急/慢性阑尾炎

十二指肠炎

脂肪肝

肠道易激综合症

急/慢性胰腺炎

//

Elabscience[®]

Elabscience Biotechnology Co.,Ltd

消化系统疾病


消化系统疾病与全身性疾病关系密切。一方面，消化系统疾病可有消化道外其它系统或全身表现，甚至在某个时期内会掩盖本系统的基本症状；另一方面，全身疾病常以消化系统症状为其主要表现或者消化道病变仅是全身疾病的一个组成部分。

近二、三十年来，由于细胞生物学、分子生物学、生物化学、内分泌学、免疫学、酶学等的快速进展，以及许多新尖精细技术的建立，使消化系统疾病的病理、病因、发病机理、诊断、防治等方面取得了很大的进展，现在的消化病学已是内科学中一门日臻完善的分支学科。

用于消化系统研究的 ELISA 试剂盒

Elabscience 根据往期文献将临床 20 种消化系统疾病（慢性胃炎 / 肠炎、药物性肝损伤、肝硬化、胰腺炎等）研究中的相关指标进行了整理，便于消化系统疾病研究人员能快速查找并获取指标信息。请访问 Elabscience 官网（www.elabscience.cn），以获取详细产品信息。

客户发表文献展示(消化系统相关)



1. Jiang S H, Li J, Dong F Y, et al. Increased Serotonin Signaling Contributes To The Warburg Effect In Pancreatic Tumor Cells Under Metabolic Stress And Promotes Growth Of Pancreatic Tumors In Mice [J]. *Gastroenterology*, 2017, 153: 277-291.
2. Raparelli V, Basili S, Carnevale R, et al. Low-Grade Endotoxemia And Platelet Activation In Cirrhosis [J]. *Hepatology*, 2017, 65(2): 571-581.
3. Yamada D, Rizvi S, Razumilava N, et al. Il-33 Facilitates Oncogene-Induced Cholangiocarcinoma In Mice By An Interleukin-6-Sensitive Mechanism[J]. *Hepatology*, 2015, 61(5): 1627-1642.
4. Mnich E, Kowalewicz-Kulbat M, Sicińska P, et al. Impact Of Helicobacter Pylori On The Healing Process Of The Gastric Barrier[J]. *World Journal of Gastroenterology*, 2016, 22(33): 7536-7558.
5. Zhan Y, Zhang F, Shi Y, et al. Potential Effects Of Calcium Binding Protein S100A12 On Severity Evaluation And Curative Effect Of Severe Acute Pancreatitis[J]. *Inflammation*, 2015, 38(1): 290-297.

消化系统疾病研究中的热点指标

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0002c	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0536c	Human AIRE(Autoimmune Regulator) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2181c	Human ALB(Albumin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2327c	Human CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H1766c	Human CFTR(Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2072c	Human CK-18/KRT18(Cytokeratin 18) ELISA Kit	6.25--400 mIU/mL	3.75 mIU/mL	Sandwich
E-EL-H0043c	Human CRP(C-Reactive Protein) ELISA Kit	0.39--25 ng/mL	0.23 ng/mL	Sandwich
E-EL-H2069c	Human CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0666c	Human CTNN β 1(Catenin, Beta 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2457c	Human CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0067c	Human FAS/CD95(Factor Related Apoptosis) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H2060c	Human GT(Gastrin) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Competitive
E-EL-H0108c	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0089c	Human IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0149c	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0099c	Human IL-2(Interleukin 2) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0102c	Human IL-6(Interleukin 6) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0103c	Human IL-10(Interleukin 10) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H2239c	Human JAK2(Janus Kinase 2) ELISA Kit	125--8000 pg/mL	75 pg/mL	Sandwich
E-EL-H2182c	Human KRAS(V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1451c	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H1946c	Human NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1359c	Human PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1970c	Human PR3(Proteinase 3) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1792c	Human PT(Prothrombin) ELISA Kit	9.38--600 ng/mL	5.63 ng/mL	Sandwich
E-EL-H1846c	Human PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1289c	Human S100A8(S100 Calcium Binding Protein A8) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H0035c	Human sCD40L(Soluble Cluster of Differentiation 40 Ligand) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0110c	Human TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H1539c	Human TLR4(Toll-Like Receptor 4) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0109c	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0910c	Human TP53(Tumor Protein p53) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H2168c	Human VWF(Von Willebrand Factor) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0965c	Human α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich

其它消化系统疾病研究中的常用指标

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0072c	Human AAP(Alanine Aminopeptidase) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0262c	Human ADA(Adenosine Deaminase) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0004c	Human ADP/Acrp30(Adiponectin) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0204c	Human AIMP1(Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1) ELISA Kit	125--8000 pg/mL	75 pg/mL	Sandwich
E-EL-H2598c	Human ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0470c	Human ApoE(Apolipoprotein E) ELISA Kit	23.44--1500 ng/mL	14.06 ng/mL	Sandwich
E-EL-H0010c	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0589c	Human BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1614c	Human BRAF(B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0017c	Human CASP3(Caspase 3) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0721c	Human CCKAR(Cholecystokinin A Receptor) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2092c	Human CCND1(Cyclin D1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0752c	Human CD4(Cluster of Differentiation 4) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H5610c	Human CD19(Cluster of Differentiation 19) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0672c	Human CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0727c	Human CENPF(Centromere Protein F) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0739c	Human CHGA(Chromogranin A) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0778c	Human COL3 α 1(Collagen Type III Alpha 1) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H1977c	Human CPT1A(Carnitine Palmitoyltransferase 1A, Liver) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H1225c	Human D2R/DRD2(Dopamine Receptor D2) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1244c	Human DLAT(Dihydrolipoyl Transacetylase) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H2485c	Human DMBT1(Deleted In Malignant Brain Tumors 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0014c	Human E-Cad(E-Cadherin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1165c	Human ELA4(Elastase 4) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich
E-EL-H1880c	Human Ep-CAM/CD326(Epithelial Cell Adhesion Molecule) ELISA Kit	1.56--100 IU/mL	0.94 IU/mL	Sandwich
E-EL-H1783c	Human ETRB(Endothelin Receptor B) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0764c	Human F5(Coagulation Factor V) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0069c	Human FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	6.25--400 ng/mL	3.75 ng/mL	Sandwich
E-EL-H0285c	Human FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	0.39--25 ng/mL	0.23 ng/mL	Sandwich
E-EL-H0068c	Human FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0375c	Human FASN(Fatty Acid Synthase) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2266c	Human FBN1(Fibrillin 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2095c	Human FLN α (Filamin A, Alpha) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0075c	Human Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0695c	Human FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2355c	Human FPN(Ferroportin) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2237c	Human GC(Glucagon) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Competitive
E-EL-H1919c	Human GHRL(Ghrelin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Competitive
E-EL-H0972c	Human GLα(Galactosidase Alpha) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H1334c	Human GSTθ1/GSTτ1(Glutathione S Transferase Theta 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0084c	Human HGF(Hepatocyte Growth Factor) ELISA Kit	125--8000 pg/mL	75 pg/mL	Sandwich
E-EL-H1659c	Human HLA-A(Leukocyte Antigen A) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H5315c	Human HSD17β10(17-Beta-Hydroxysteroid Dehydrogenase Type 10) ELISA Kit	25--1600 pg/mL	15 pg/mL	Sandwich
E-EL-H0159c	Human IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2532c	Human IFN-α(Interferon Alpha) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0085c	Human IFN-β(Interferon Beta) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0088c	Human IL-1α(Interleukin 1 Alpha) ELISA Kit	3.91--250 pg/mL	2.34 pg/mL	Sandwich
E-EL-H0107c	Human IL-23(Interleukin 23) ELISA Kit	39.06--2500 pg/mL	23.44 pg/mL	Sandwich
E-EL-H1030c	Human IL-2Rα(Interleukin 2 Receptor Alpha) ELISA Kit	93.75--6000 pg/mL	56.25 pg/mL	Sandwich
E-EL-H0101c	Human IL-4(Interleukin 4) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0191c	Human IL-5(Interleukin 5) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0424c	Human IL-10Rβ(Interleukin 10 Receptor Beta) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H1545c	Human IRF5(Interferon Regulatory Factor 5) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2506c	Human ITGα4(Integrin Alpha 4) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1895c	Human L1CAM/CD171(L1-Cell Adhesion Molecule) ELISA Kit	93.75--6000 pg/mL	56.25 pg/mL	Sandwich
E-EL-H2456c	Human LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0113c	Human LEP(Leptin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H5200c	Human LTF/LF(Lactoferrin) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H5317c	Human M I CA(major histocompatibility complex class I chain-related gene A) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1770c	Human MAdCAM1(Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0020c	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H1445c	Human MMP-2(Matrix Metalloproteinase 2) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0844c	Human MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H2055c	Human MTL(Motilin) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Competitive
E-EL-H0857c	Human MYC(V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1894c	Human NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H2104c	Human PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2214c	Human PDGFRα(Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2245c	Human PL(Pancreatic Lipase) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H1360c	Human PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1196c	Human PRSS1(Protease, Serine, 1) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1986c	Human RAP1A(Ras-related protein Rap-1A) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0355c	Human RARS(Arginyl tRNA Synthetase) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H2243c	Human REG4(Regenerating Islet Derived Protein 4) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1379c	Human RNASE3/ECP(Ribonuclease A3/Eosinophil Cationic Protein) ELISA Kit	39.06--2500 pg/mL	23.44 pg/mL	Sandwich
E-EL-H1290c	Human S100A9(S100 Calcium Binding Protein A9) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1151c	Human SCT(Secretin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Competitive
E-EL-H0930c	Human SHH(Hedgehog Homolog, Sonic) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0135c	Human sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0933c	Human Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H5540c	Human SPINK1(Serine Peptidase Inhibitor Kazal Type 1) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H1923c	Human SST(Somatostatin) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Competitive
E-EL-H2130c	Human STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2132c	Human STAT3(Signal Transducer and Activator of Transcription 3) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0171c	Human TF(Transferrin) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H1657c	Human TGF β 1(Transforming Growth Factor Beta Induced Protein) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H0184c	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1588c	Human TPO(Thrombopoietin) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0439c	Human TTN(Titin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2342c	Human TTR(Transthyretin) ELISA Kit	15.63--1000 ng/mL	9.38 ng/mL	Sandwich
E-EL-H1795c	Human VCL(Vinculin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H2155c	Human VIP(Vasoactive Intestinal Peptide) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Competitive
E-EL-H0070c	Human α FP(Alpha-Fetoprotein) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0980c	Human α -Glu(Alpha Glucosidase) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1012c	Human γ GT1(Gamma Glutamyltransferase 1) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich



网 址: www.elabscience.cn

电 话: 027-87526315

E-mail: techsupport@elabscience.cn tech@elabscience.cn

地 址: 武汉市东湖新技术开发区高新大道 858 号生物医药园二期 B18 栋

