


ELISA Kits

for Cardiovascular Medicine Research



//

病态窦房节综合症
风湿性心脏病
高血压
冠心病
结节性多动脉炎
雷诺氏病
心房纤颤
心肌病
心肌梗塞
心力衰竭
心律失常

//

Elabscience®

Elabscience Biotechnology Co.,Ltd

心血管内科疾病

心内科，即心血管内科，是各级医院大内科为了诊疗心血管疾病而设置的一个临床科室，治疗的疾病包括心绞痛、高血压、猝死、心律失常、心力衰竭、早搏、心律不齐、心肌梗死、心肌病、心肌炎、心肌梗塞等心血管疾病。

心血管疾病是全球最常见的死因之一，除了非洲之外心血管疾病在死因排行中都名列前茅，成为严重威胁人民健康的主要疾病，所有心血管疾病的成因、防止及治疗都是医学研究的重要范畴，有超过数百个的科学研究在每星期刊登。

用于心血管内科疾病研究的 ELISA 试剂盒

Elabscience 根据往期文献将临床 11 种心血管疾病（高血压、冠心病、心肌梗塞、风湿性心脏病等）涉及指标进行了整理，便于心血管内科疾病研究人员能快速查找并获取指标信息。请访问 Elabscience 官网（www.elabscience.cn），以获取详细产品信息。

客户发表文献展示（心血管系统相关）

1. Huang B, Wang W H, Li Q C, et al. Osteoblasts Secrete Cxcl9 To Regulate Angiogenesis In Bone [J]. Nature Communications, 2016.
2. Liu Y, Zhong G, Zhao Z. Effect Of Chronic Left Ventricular Unloading On Myocardial Remodeling: Multimodal Assessment Of Two Heterotopic Heart Transplantation Techniques [J]. Journal of Heart and Lung Transplantation, 2014,34(4): 594-603.
3. Sun L, Yang X X, Li Q, et al. Activation Of Adiponectin Receptor Regulates Proprotein Convertase Subtilisin/Kexin Type 9 Expression And Inhibits Lesions In ApoE-Deficient Mice [J]. Arteriosclerosis Thrombosis and Vascular Biology, 2017, 37: 1290-1300.
4. Wu M, Claus P, Driessche N V, et al. Placental Growth Factor 2-A Potential Therapeutic Strategy For Chronic Myocardial Ischemia [J]. International Journal of Cardiology, 2016, 203: 534-542.
5. He W, Ye S, Zeng C, et al. Hypothermic Oxygenated Perfusion (HOPE) Attenuates Ischemia/Reperfusion Injury In The Liver Through Inhibition Of The TXNIP/NLRP3 Inflammasome Pathway In A Rat Model Of Donation After Cardiac Death [J]. The FASEB Journal, 2018: fj. 201800028RR.




心血管内科疾病研究中的热点指标

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0002c	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0257c	Human ACTN2(Actinin Alpha 2) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0300c	Human AGT(Angiotensinogen) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H2181c	Human ALB(Albumin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0329c	Human ANG II R1(Angiotensin II Receptor 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0532c	Human ANP(Atrial Natriuretic Peptide) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Competitive
E-EL-H0125c	Human ApoA1(Apolipoprotein A1) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich
E-EL-H0470c	Human ApoE(Apolipoprotein E) ELISA Kit	23.44--1500 ng/mL	14.06 ng/mL	Sandwich
E-EL-H0598c	Human BNP(Brain Natriuretic Peptide) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Competitive
E-EL-H0043c	Human CRP(C-Reactive Protein) ELISA Kit	0.39--25 ng/mL	0.23 ng/mL	Sandwich
E-EL-H0649c	Human cTn-I/TNNT3(cardiac Troponin I) ELISA Kit	0.39--25 ng/mL	0.23 ng/mL	Sandwich
E-EL-H0646c	Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H1505c	Human Des(Desmin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1850c	Human DSP(Desmoplakin) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0064c	Human ET-1(Endothelin 1) ELISA Kit	1.25--80 pg/mL	0.75 pg/mL	Sandwich
E-EL-H0764c	Human F5(Coagulation Factor V) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H2266c	Human FBN1(Fibrillin 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0882c	Human GATA4(GATA Binding Protein 4) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0108c	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0149c	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0102c	Human IL-6(Interleukin 6) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0103c	Human IL-10(Interleukin 10) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0844c	Human MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0717c	Human MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0753c	Human NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0755c	Human NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H1947c	Human PLN(Phospholamban) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0119c	Human REN(Renin) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0035c	Human sCD40L(Soluble Cluster of Differentiation 40 Ligand) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0110c	Human TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0166c	Human TM(Thrombomodulin) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H0109c	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0439c	Human TTN(Titin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2168c	Human VWF(Von Willebrand Factor) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich

其它心血管内科疾病研究中的常用指标

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0464c	Human ApoB(Apolipoprotein B) ELISA Kit	78.13--5000 ng/mL	46.88 ng/mL	Sandwich
E-EL-H0472c	Human ApoH(Apolipoprotein H) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0418c	Human BAG3(Bcl2 Associated Athanogene 3) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0010c	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0590c	Human BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1614c	Human BRAF(B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H5495c	Human CECR1(Cat Eye Syndrome Chromosome Region, Candidate 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1519c	Human CRYαB(Crystallin Alpha B) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1436c	Human DMD(Dystrophin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2390c	Human DPP6(Dipeptidyl Peptidase VI) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich
E-EL-H0062c	Human ENG(Endoglin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H3641c	Human F8(Coagulation Factor VIII) ELISA Kit	46.88--3000 mIU/mL	28.13 mIU/mL	Sandwich
E-EL-H0074c	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0883c	Human GATA5(GATA Binding Protein 5) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1888c	Human GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1925c	Human Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0085c	Human IFN-β(Interferon Beta) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H2665c	Human INS(Insulin) ELISA Kit	0.78--50 μIU/mL	0.47 μIU/mL	Sandwich
E-EL-H2203c	Human ITGB3(Integrin Beta 3) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0113c	Human LEP(Leptin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H1446c	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1451c	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0142c	Human MYO(Myoglobin) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich
E-EL-H1946c	Human NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1368c	Human NUP155(Nucleoporin 155kDa) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2104c	Human PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2205c	Human PKP2(Plakophilin 2) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1014c	Human POMC(Pro-Opiomelanocortin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2298c	Human PON1(Paraoxonase 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1970c	Human PR3(Proteinase 3) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1792c	Human PT(Prothrombin) ELISA Kit	9.38--600 ng/mL	5.63 ng/mL	Sandwich
E-EL-H1843c	Human PTGIS(Prostaglandin I Synthase) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0940c	Human SAG(S Antigen) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0135c	Human sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H1113c	Human SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Competitive
E-EL-H0876c	Human sSELE(Soluble E-Selectin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H2459c	Human TNNC1(Troponin C Type 1) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H1954c	Human TNNI3K(Serine/threonine-protein kinase TNNI3K) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H2106c	Human tPA(Tissue-type Plasminogen Activator) ELISA Kit	0.63--40 ng/mL	0.38 ng/mL	Sandwich
E-EL-H2342c	Human TTR(Transthyretin) ELISA Kit	15.63--1000 ng/mL	9.38 ng/mL	Sandwich
E-EL-H1795c	Human VCL(Vinculin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0980c	Human α -Glu(Alpha Glucosidase) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich



网 址: www.elabscience.cn

电 话: 027-87526315

E-mail: techsupport@elabscience.cn tech@elabscience.cn

地 址: 武汉市东湖新技术开发区高新大道 858 号生物医药园二期 B18 栋

